



AutoClean[®]

Washer Disinfector





AutoClean® latest generation washer-disinfector.

Its 19L tank allows it to manage a large instrumentation capacity in one go.

The **AutoClean** allows you to have a use adapted to your needs thanks to its half-load function.

Useful capacity 15L at full load and 9L at half load.

Fully automatic, it saves valuable time for soaking, cleaning and drying work.

Thanks to its design and modern design, the AutoClean will fit into any dental practice.

Its touch control panel is intuitive, it allows consultation of the history of use and also gives the possibility of manual use of each sequence.

- Automatic or manual cycle
- Economic use (product concentration 1%)
- Optimize your assistant's time
- Reduce the congestion


Intuitive control panel



Activation of full load

Activation of half load

Visualization position in cycle



Technical characteristics 19L

- Supply: 230 V AC - 50Hz
- Ultrasonic power: 400W
- Number of transducers: 8
- Cleaning frequency: 40 KHz
- Drying power: 1000W
- Heating power: 200W
- Dimensions of the tank (mm): 325 X 300 X 200 mm
- Dimensions on table (L x l x h): 520 X 490 X 80 mm
- Dimensions under the table (L x l x h): 420 X 430 X 340 mm
- Dimensions of the cut (L x l x h): 440 x 450 mm
- USB: Historical cycles report
- Drain: Motorized valve
- Function: Manual / Automatic
- Lid lock: Stop cycle
- Volume of the tank: 19L
- Useful volume full load: 15L
- Useful volume half load: 9L
- Color touch screen: 10 cm
- Weight: 35 kg
- Delivered with basket

Reference	Designation
MHC190A	Automatic Cleaning machine

Technical characteristics 13L

- Supply: 230 V AC - 50Hz
- Ultrasonic power: 300W
- Number of transducers: 6
- Cleaning frequency: 40 KHz
- Drying power: 1000W
- Heating power: 200W
- Dimensions of the tank (mm): 325 X 300 X 150 mm
- Dimensions on table (L x l x h): 520 X 490 X 80 mm
- Dimensions under the table (L x l x h): 420 X 430 X 340 mm
- Dimensions of the cut (L x l x h): 440 x 450 mm
- USB: Historical cycles report
- Drain: Motorized valve
- Function: Manual / Automatic
- Lid lock: Stop cycle
- Volume of the tank: 13L
- Useful volume full load: 9L
- Useful volume half load: 6L
- Color touch screen: 10 cm
- Weight: 32k g
- Delivered with basket

Reference	Designation
MHC130A	Automatic Cleaning machine





AutoClean NET

Instrumentation disinfectant detergent.
Cleaning and disinfection of medical and surgical instruments and equipment.
Cleaning in automatic ultrasonic tanks.
5L Container

Reference	Designation
MHC190B	Decontamination solution AutoClean NET



AutoClean RN

Neutralization of alkaline residues from detergents and water during the automatic washing cycle of medical-surgical instruments and medical devices.
5L Container

Reference	Designation
MHC190C	Cleaning solution AutoClean RN



Half basket for AutoClean automatic cleaner

Reference	Designation
MHC190F	Half basket / piece



Basket for AutoClean automatic cleaner

Reference	Designation
MHC190G	Basket AutoClean





MHC Technology

1001 chemin Hugues Bérenguier

06610 LA GAUDE (France)

Tel: 04.93.58.17.16. Fax: 04.93.58.14.75

email: info@mhc-technology.com

www.mhc-technology.com